BonOs® R NF Genta

Bone Cement

Instructions for use	2	en
Gebrauchsanweisung	7	de
Notice d'instruction	12	fr
Instrucciones de utilización	17	es
Istruzioni per l'uso	22	it
Instruções de utilização	27	pt
Návod k použití	32	cs
Gebruiksaanwijzing	36	nl
تعليمات للإستخدام	41	ar
사용 설명서	47	ko
Инструкция по применению	51	ru
Kullanım talimatları	56	tr

checked QA
26. Nov. 2018
sign: *N. S.*

OSARTIS GmbH
Lagerstraße 11-15
64807 Dieburg
Germany
phone: +49 (0) 6071 – 929 0
fax: +49 (0) 6071 – 929 100
e-mail: info@osartis.de
web: www.osartis.de

CE 0123

File: 135-0576-03

Purpose and properties

BonOs® R NF Genta is a fast-setting acrylic resin with addition of gentamicin sulfate for use in bone surgery. Mixing the two separate sterile components produces a ductile bone cement which, after hardening, fixes the implant and transfers stresses produced during movement evenly to the bone. The added antibiotic, gentamicin sulfate, protects the implant and the surrounding tissue from colonization with pathogens that are sensitive to gentamicin. BonOs® R NF Genta cement powder also contains insoluble zirconium dioxide as an X-ray contrast medium. BonOs® R NF Genta does not emit a signal and does not pose a safety risk in a magnetic resonance environment.

Composition BonOs® R NF Genta

The cement powder contains	%	1 x 20	1 x 40	1 x 60
Poly (methylacrylate/ Methyl methacrylate)	82.48	16.82 g	33.65 g	50.47 g
Zirconium dioxide	14.70	3.00 g	6.00 g	9.00 g
Benzoyl peroxide	0.86	0.18 g	0.35 g	0.53 g
Gentamicin sulfate (Gentamicin base)	1.96 (1.23)	0.40 g (0.25 g)	0.80 g (0.50 g)	1.20 g (0.75 g)

The liquid component contains	%	1 x 20	1 x 40	1 x 60
Methyl methacrylate (stabilized with 60 ppm hydroquinone)	98.47	9.26 g	18.51 g	27.77 g
N,N-dimethyl-p-toluidine	1.53	0.14	0.29 g	0.43 g

BonOs® R NF Genta is offered in the following pack sizes

Pack size	Powder mass	Liquid volume
2 x 20	2 x 20.4 g	2 x 10 ml
1 x 40	1 x 40.8 g	1 x 20 ml
2 x 40	2 x 40.8 g	2 x 20 ml
1 x 60	1 x 61.2 g	1 x 30 ml

Indications

BonOs® R NF Genta is indicated for fixation of synthetic resin and metal prosthesis components during partial or total replacement of the hip and knee as well as other joints, if an infection with organisms that are sensitive to gentamicin is diagnosed or suspected. The antibiotic provides protection against bacterial colonization of the implant and the surrounding tissue.

Contraindications

The use of BonOs® R NF Genta is contraindicated, if muscle wasting or neuromuscular compromise in the affected limb renders the procedure unjustifiable. BonOs® R NF Genta should not be used in the event of known hypersensitivity to any of its constituents or in patients with severe renal failure.

Information for use

Before using BonOs® R NF Genta for the first time, surgeons should familiarize themselves with the mixing and application process. A trial run of the mixing process is recommended. If special mixing and application techniques are used, the surgeon must read the relevant instructions first.

The protective outer packaging (aluminum/PE pouch) and blister pack containing the ampoule should be removed from the carton in the non-sterile area. After removal from the protective outer packaging, the polyethylene pouch ("peel-off pouch"), the outside of which is not sterile, and the blister pack containing the ampoule are opened one after the other under strict aseptic technique by a circulating nurse and passed to a member of the surgical team in the sterile area. After careful preparation of the implant site the cement can be applied to the bone, using a cement syringe or with other application techniques (details of how to do this can be found in the instructions for the system used). The ampoule is opened by breaking its neck and the inner pouch containing the cement powder cut open using sterile scissors.

Dosage

The amount of BonOs® R NF Genta required depends on the anatomical conditions of the patient and the implant used. If a large amount is needed, additional packs (powder and monomer) may be mixed together. However, at least one complete unit (the contents of one pouch and one ampoule) must always be mixed. In total, no more than 160 g should be applied at once. It is advisable to keep additional packs of BonOs® R NF Genta ready as a precaution.

Tools required for mixing BonOs® R NF Genta

Sterile working area, sterile porcelain bowls, stainless steel bowls or plastic bowls that are suitable for monomers, sterile mixing spoons or spatulas made of porcelain or high-grade steel, or a sterile mixing system for bone cements. To utilize a modern cementing technique, the use a vacuum mixing system is recommended.

Manual mixing and application

To make the mixture, empty the entire contents of the required ampoules into a suitable sterile, inert mixing vessel. Then add the

entire contents of the corresponding number of packs of powder to the liquid and, using a suitable spatula, carefully mix the components for about 30 seconds until the powder has completely absorbed the liquid and a homogenous dough has formed. The cement can be applied when it stops sticking to the surgeon's gloves and has reached the desired viscosity. The application phase ends when the dough becomes rubbery and elastic and no longer binds completely when kneaded. If application of the cement continues, even filling of the bone cannot be ensured and there is a risk that the implant will loosen prematurely (see Fig. 1).

Vacuum mixing

In order to reduce porosity, the cement can be mixed in a vacuum mixing system. BonOs® R NF Genta is designed for handling between 17°C and 25°C, however, prechilling is recommended. The mixing time should be approximately 30 seconds. The contents of the required ampoules are emptied into the mixing vessel and the corresponding number of powder units added. The further procedure can be found in the instructions, supplied by the manufacturer of the mixing system (see Fig. 2).

Prechilling

Prechilling is recommended if a lower viscosity or a prolonged handling time of the bone cement is required. Prechilling of the bone cement for at least 24 hours makes mixing more convenient and reduces the viscosity of the bone cement. The mixing time is also 30 seconds, however the application and hardening phase is longer. With decrease in temperature, viscosity decreases and handling as well as hardening time increase. For handling characteristics under vacuum mixing at for example 4 °C please refer to fig. 3. Do not prechill BonOs® R NF Genta at temperatures below 4°C.

Application with a cement syringe

For several minutes after mixing, the cement can be applied with a syringe, although the process must be monitored carefully by the surgeon due to the increasing viscosity over time. In hip joint replacement, the use of a restrictor or plug in the femoral canal is strongly recommended.

Notes on use

- For a good fixation the implant should be inserted during the application phase and held in place until the cement hardens.
 - Surplus cement should be removed before it hardens.
 - The temperature-time graphs should be noted.
 - Prechilling the cement components reduces viscosity and prolongs handling and hardening time.
 - Prechilling of the cement components at 4°C is recommended when mixed with a vacuum mixing device.
 - Handling time and polymerization are heavily dependent on the temperature of the components and environment. Hardening time is reduced by higher temperatures and prolonged by lower temperatures. Viscosity increases with the progression of polymerization, i.e. the duration of the handling phase.
 - The addition of any other powders or liquids can reduce solidity and/or impair handling characteristics and must therefore be avoided.
 - During mixing and application, it is important to minimize air entrapment.
 - The powder and liquid components have been carefully designed to complement each other. The entire contents of the pouch and ampoule must always be mixed.
- USING ONLY PARTS OF THE COMPONENTS IS NOT PERMITTED!**
- It is recommended to verify the correct implantation with suitable imaging procedures.
 - If the cement is applied in a too low viscous state, the bleeding pressure can intrude into the cement mass and can lower the longevity of the implant.

Use in joint surgery

When using BonOs® R NF Genta in joint surgery, a modern cementing technique should be used to limit undesirable effects. A prerequisite for this is careful preparation of the implant site with thorough rinsing (e.g. pulsatile lavage) and drying before application of the cement. In order to prevent pressure from building in the medullary canal during implantation, adequate drainage is recommended. Further prerequisites for better fixation of the implant include filling the entire medullary canal with cement using a restrictor, completely surrounding the implant with a cement mantle (ideally 2-5 mm thick) and achieving the ideal biomechanical fit in the bone.

Adverse events

After preparation of the implant site and immediately after application of cement and implantation, the rise in pressure in the medullary canal can lead to a temporary drop in blood pressure. In rare cases, pulmonary embolism and myocardial infarction are also observed. These cardiovascular and respiratory side effects, which are known as implantation syndrome, result mainly from infiltration of bone marrow constituents into the venous system.

The following additional undesirable effects have occurred after using acrylate cements: temporary reduction in blood pressure, elevated serum levels of gamma-glutamyl transferase (gamma-GT) up to 10 days after surgery, thrombophlebitis, hemorrhage and hematoma, loosening or dislocation of the implant, superficial or deep wound infection, trochanteric bursitis, heterotopic ossification and trochanteric detachment, cardiovascular reactions such as temporary heart rhythm disorders, short-term cardiac conduction irregularities, arrhythmia, myocardial infarction and cardiac arrest, hypoxemia, bronchospasm, pulmonary embolism, apoplexy. There have been rare reports of hypotension with anaphylaxis including anaphylactic shock associated with cardiac arrest and sudden death.

Other adverse events that may be attributable to the use of PMMA bone cements are: allergic pyrexia, hematuria, dysuria, bladder fistulas, local neuropathy and vascular erosion or occlusion as well as postoperative irritation of the sciatic nerve, because bone cement has been placed outside the intended area of application.

In isolated cases, the gentamicin contained in BonOs® R NF Genta can cause hypersensitivity reactions. In principle, the typical side effects of using gentamicin, particularly hearing problems and kidney damage, cannot be ruled out completely.

However, these side effects are extremely unlikely to occur, because of the very low serum level of gentamicin (<1 µg/ml). Monomer vapors can irritate the respiratory tract and eyes and may damage organs.

Interactions

The tendency of gentamicin to block neuromuscular transmission can be intensified by concomitant administration of muscle relaxants, e.g. D-tubocurarine, suxamethonium or pancuronium, as well as by ether. Concomitant administration of potentially neurotoxic and/or nephrotoxic substances, e.g. cisplatin, other aminoglycosides, streptomycin, cefaloridine, viomycin, polymyxin B or polymyxin E, can increase the toxicity of gentamicin. However, interaction is very unlikely to occur because of the low serum levels of gentamicin.

Warnings and precautions

Use during pregnancy and lactation

There are no adequate studies of the use of acrylate cements during pregnancy and lactation or of their effect on human fertility. During pregnancy and lactation, the surgeon should weigh the benefit for the mother against the potential risk to the child before using BonOs® R NF Genta bone cement.

Warnings relating to the age of the patients to be treated

There are no adequate studies on the use of acrylic bone cements in children. As the possibility of acrylate cements having adverse effects on bone growth cannot be ruled out, the use of BonOs® R NF Genta in children and patients who are still growing is inadvisable.

As reported in the clinical literature, in younger patients, cementless components perform as well or better than cemented components.

During and immediately after application of the bone cement/implantation, blood pressure, pulse and respiration must be monitored carefully, and appropriate measures taken, if significant changes occur. If the patient develops pulmonary or cardiovascular symptoms, appropriate monitoring of blood loss is required. In the event of acute respiratory failure, anesthesiological measures should be initiated immediately.

Hypotensive reactions have occurred between 10 and 165 seconds following application to bone; they have lasted from 30 seconds to 5 or more minutes. Some have progressed to cardiac arrest.

Polymerization (hardening) of BonOs® R NF Genta is an exothermic reaction. The heat given off during this reaction may damage bone or other tissue in the region of the implant.

The training and experience of the surgeon are very important when handling acrylate bone cements. Instructions for handling and mixing the cement and preparing the implantation site must be followed carefully. Before using BonOs® R NF Genta, the surgeon must be thoroughly familiar with its properties and handling characteristics. As the handling and setting characteristics of BonOs® R NF Genta depend on the temperature and mixing technique, they are best determined by the surgeon on the basis of actual experience. For this reason, the surgeon is strongly recommended to carry out a trial run of the entire mixing, handling and setting process before performing a surgical procedure with BonOs® R NF Genta. Inadequate fixation and unexpected postoperative events impair the interface between the cement and bone. This can cause micro-movements, which may lead to the formation of a layer of fibrous tissue and premature implant failure. Premature loosening of the implant in the cement mantle is also possible. Therefore, long-term regular follow-up examinations are recommended for all patients. Whenever clinically indicated additional appropriate systemic antibiotic coverage should be used on the day of surgery. Careful consideration should be given to the potential additional safety aspects of the additional antibiotics.

Methyl methacrylate is a volatile, flammable liquid. The vapors produced during the mixing process can irritate the respiratory tract and eyes and cause general malaise and headache. Such symptoms can be reduced with adequate ventilation or by using closed mixing systems. The monomer (methyl methacrylate) is lipid-soluble. Direct skin contact with the liquid monomer should be avoided as far as possible, as allergic reactions (contact dermatitis) cannot be ruled out. It is therefore advisable to wear an additional pair of polyethylene (PE) gloves under the normal surgical gloves when handling the cement. The following materials have also proved to be suitable for protective gloves: PVP (polyethylene, ethylene vinyl alcohol, polyethylene) and viton butyl.

Contact lenses should be protected from the escaping monomer vapors.

The time point when the prosthesis can be fully loaded, as well as precautions and actions to be avoided after surgery depends on the surgical procedure, the type of associated implant and the condition of the patient and should be judged by the responsible surgeon. Medical staff must inform the patient of all precautions to be taken.

Removal of bone cement in case of revision

For revision, a radical debridement that entitles removal of cement, and potentially infected and devitalized tissue is performed. Intramedullary reaming of the canal followed with copious use of saline pulse lavage is recommended. The procedure to remove the implant is a case by case decision by the responsible surgeon depending on the condition of the implant and the patient.

Storage

BonOs® R NF Genta should be stored away from direct sunlight. BonOs® R NF Genta is stored at temperatures between 0°C (32°F) and +25°C (77°F). Do not use BonOs® R NF Genta after the expiry date.

Shelf life/sterility

The expiry date is printed on the batch label of the carton box, on the patient label, the aluminum pouch and inner pouch of the cement powder and the label of the monomer blister lid as well as the monomer ampoule. BonOs® R NF Genta must not be used

after this date. The contents of unused open or damaged packs must not be resterilized and should therefore be discarded. The polymer powder must not be used, if it shows yellow discoloration. The liquid monomer (aseptically filled), the ampoule itself (outside sterilized using ethylene oxide), the inside of the blister (sterilized using ethylene oxide), the cement powder and the inner PE/ paper pouch are all sterile. The product must not be used, if the packaging is damaged.

Single-use

BonOs® R NF Genta must never be reused. Due to the functionality (hardening) of PMMA bone cements, BonOs® R NF Genta is only suitable for use within the specified application period. One unit of BonOs® R NF Genta must be used for only one single patient.

Safe disposal

1. Mixed cement should set before it is discarded with hospital waste.
2. For separate disposal of powder and liquid, please ask local waste disposal authorities.

Information



For further information, please contact your supplier or the manufacturer directly.

Symbol definitions

	Manufacturer		Consult instructions for use
	Use-by date		Caution
	Batch code		CE marking of conformity
	Catalogue number		Registered trademark
	Sterilized using aseptic processing techniques		Flammable
	Sterilized using ethylene oxide		Harmful
	Do not re-sterilize		Health Hazard
	Do not use if package is damaged		Ampoule
	Keep away from sunlight		Blister
	Temperature limit		Pouch
	Do not re-use		


Mixing phase: A Waiting phase: B

Fig. 1: Manual mixing

Fig. 3: Prechilled vacuum mixing¹

¹ Different mixers may have different effects on the properties of the bone cement

Application phase: C Hardening phase: D

Fig. 2: Vacuum mixing¹

OSARTIS GmbH
Lagerstraße 11-15
64807 Dieburg
Germany
phone: +49 (0) 6071 929 0
fax: +49 (0) 6071 929 100
e-mail: info@osartis.de
www.osartis.de

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