BonOs® R NF

Bone Cement

checked QA

09. Dez. 2019

sign:

OSARTIS GmbH

Lagerstraße 11-15

64807 Dieburg

Germany

phone: +49 (0) 6071 – 929 0

fax: +49 (0) 6071 – 929 100

e-mail: info@osartis.de

web: www.osartis.de

File: 135-0575-03

Instructions for use	2	en
Gebrauchsanweisung	6	de
Notice d'instruction	10	fr
Instrucciones de utilización	14	es
Istruzioni per l'uso	18	it
Instruções de utilização	22	pt
Návod k použití	26	cs
Brugsanvisning	30	da
Οδηγίες χρήσης	34	el
Kasutusjuhend	38	et
Käyttöohje	42	fi
Upute za uporabu	46	hr
Használati utasítás	50	hu
Naudojimo instrukcija	54	lt
Lietošanas pamācība	58	lv
Gebruiksaanwijzing	62	nl
Instrukcje użytkowania	66	pl
Návod na použitie	70	sk
Bruksanvisning	74	sv
مادخت سالات ادادشرا	78	ar
사용 설명서	82	ko
Инструкция по применению	86	ru
Kullanım talimatları	90	tr

Purpose and properties

BonOs® R NF is a fast-setting acrylic resin for use in bone surgery. Mixing the two separate sterile components produces a ductile bone cement which, after hardening, fixes the implant and transfers stresses produced during movement evenly to the bone. BonOs® R NF cement powder also contains insoluble zirconium dioxide as an X-ray contrast medium. BonOs® R NF does not emit a signal and does not pose a safety risk in a magnetic resonance environment. The light colour of BonOs® R NF makes it easy to identify.

Composition BonOs® R NF

The cement powder contains		
Poly (methacrylate/methyl methacrylate)	%	1 x 40
Zirconium dioxide	84.1 %	33.65 g
Benzoyl peroxide	15.0 %	6.00 g
	0.9 %	0.35 g
The liquid component contains		
Methyl methacrylate (stabilized with 60 ppm hydroquinone)	%	1 x 40
N,N-dimethyl-p-toluidine	98.47 %	18.51 g
	1.53 %	0.29 g

BonOs® R NF is offered in the following pack sizes

Pack size	Powder mass	Liquid volume
1 x 40	1 x 40 g	1 x 20 ml
2 x 40	2 x 40 g	2 x 20 ml

Indications

In partial or total replacement of the hip, knee and other joints, BonOs® R NF is indicated for fixation of synthetic resin and metal prosthesis components in uninfected vital bones if joint reconstruction is required.

In tumour surgery, BonOs® R NF is used in combination with internal fixation to fill bone cavities after tumour removal.

Contraindications

The use of BonOs® R NF is contraindicated if muscle wasting or neuromuscular compromise in the affected limb renders the procedure unjustifiable. BonOs® R NF should not be used in the event of known hypersensitivity to any of its constituents or of active or inadequately treated infections in the region of the implant.

Use during pregnancy and lactation

There are no adequate studies of the use of acrylic cements during pregnancy and lactation or of their effect on human fertility. During pregnancy and lactation the surgeon should weigh the benefit for the mother against the potential risk to the child before using BonOs® R NF bone cement.

Warnings relating to the age of the patients to be treated

There are no adequate studies on the use of acrylic cements in children. As the possibility of acrylic cements having adverse effects on bone growth cannot be ruled out, the use of BonOs® R NF in children and patients who are still growing is inadvisable. In younger patients (< 50 years old), if possible, joint replacement procedures that do not use cement should be favoured over those that use acrylic cements for securing the implant.

Information for use

Before using BonOs® R NF for the first time, surgeons should familiarise themselves with the mixing and application process. A trial run of the mixing process is recommended. If special mixing and application techniques are to be used, the surgeon must read the relevant instructions first.

The protective outer packaging (aluminium/PE pouch) and blister pack containing the ampoule should be removed from the carton in the non-sterile area. After removal from the protective outer packaging, the polyethylene pouch ("peel-off pouch"), the outside of which is not sterile, and the blister pack containing the ampoule are opened one after the other in sterile conditions by a circulating nurse and passed to a member of the surgical team in the sterile area. After careful preparation of the medullary canal, the cement can be applied to the bone manually, using a cement syringe or with other application systems (details of how to do this can be found in the instructions for the system used). The ampoule is opened by breaking its neck and the inner pouch containing the cement powder cut open using sterile scissors.

Dosage

The amount of BonOs® R NF required depends on the anatomical conditions of the patient and the implant used. If a large amount is needed, additional packs (powder and monomer) may be mixed together. However, at least one complete unit (the contents of one pouch and one ampoule) must always be mixed. In total, no more than 4 packs should be mixed at once. It is advisable to keep additional packs of BonOs® R NF ready as a precaution.

Tools required for mixing BonOs® R NF

Sterile working area, sterile porcelain bowls, stainless steel bowls or plastic bowls that are suitable for monomers, sterile mixing

spoons or spatulas made of porcelain or high-grade steel, or a sterile mixing system for bone cements. For state-of-the-art technology, use a vacuum mixing system.

Manual mixing and application

To make the mixture, empty the entire contents of the required ampoules into a suitable sterile, inert mixing vessel. Then add the entire contents of the corresponding number of packs of powder to the liquid and, using a suitable spatula, carefully mix the components for about 30 seconds until the powder has completely absorbed the liquid and a homogenous mixture has formed. The cement can be applied when it stops sticking to the surgeon's gloves and has reached the desired viscosity. The application phase ends when the paste becomes rubbery and elastic and no longer binds completely when kneaded. If application of the cement continues, even filling of the bone cannot be ensured and there is a risk that the implant will loosen prematurely (see Fig. 1).

Vacuum mixing

In order to reduce porosity, the cement can be mixed in a vacuum mixing system. The prechilling is not mandatorily necessary. BonOs® R NF is designed for handling between 17°C and 25°C, however, prechilling is recommended. The mixing time should be approximately 30 seconds. The contents of the required ampoules are placed in the mixing vessel and the corresponding number of powder units added. The further procedure can be found in the instructions supplied by the manufacturer of the mixing system (see Fig. 2).

Prechilling

Prechilling of the bone cement makes mixing more convenient and reduces the viscosity of the bone cement. To obtain optimal results store the cement components at 4°C for at least 24 hours. The mixing time is also 30 seconds, however the application phase is longer. Prechilling is recommended for vacuum mixing (see Fig. 3).

Application with a cement syringe

For several minutes after mixing, the cement can be applied with a syringe, although the process must be monitored carefully by the surgeon. In hip joint replacement, the use of a restrictor or plug in the femoral canal is strongly recommended.

Important general information

- For good fixation of the implant, it should be inserted during the application phase and held in place until the cement hardens.
- Surplus cement should be removed before it hardens.
- The temperature-time graphs should be noted.
- Pre-chilling the cement components reduces viscosity and prolongs handling and hardening time.
- Prechilling of the cement components at 4°C is recommended when mixed with a vacuum mixing device.
- Handling time and polymerisation are heavily dependent on the temperature of the components and environment. Hardening time is reduced by higher temperatures and prolonged by lower temperatures. Viscosity increases with the progression of polymerisation, i.e. the duration of the handling phase.
- The addition of any other powders or liquids can reduce solidity and impair handling characteristics.
- During mixing and application it is important to minimize air entrapment.
- The powder and liquid components have been carefully designed to complement each other. The entire contents of the pouch and ampoule must always be mixed.

USING ONLY PART OF THE COMPONENTS IS NOT PERMITTED!

10. It is recommended to verify the correct implantation with suitable imaging procedures.

Use in joint surgery

When using BonOs® R NF in joint surgery, a modern cementing technique should be used in order to limit undesirable effects and ensure stable long-term fixation of the implant. A prerequisite for this is careful preparation of the implant site with thorough rinsing (e.g. pulsatile lavage) and drying before application of the cement. In order to prevent pressure from building in the medullary canal during implantation, adequate drainage is recommended. Further prerequisites for better fixation of the implant include filling the entire medullary canal with cement using a restrictor, completely surrounding the implant with a cement mantle (ideally 2-5 mm thick) and achieving the ideal biomechanical fit in the bone.

Use in tumour surgery

After removal of the tumour, bone cavities are filled with BonOs® R NF. Plates or screws are then inserted to ensure the necessary stability.

Side effects

After preparation of the implant site and immediately after application of cement and implantation, the rise in pressure in the medullary canal can lead to a temporary drop in blood pressure. In rare cases, pulmonary embolism and myocardial infarction are also observed. These cardiovascular and respiratory side effects, which are known as implantation syndrome, result mainly from infiltration of bone marrow constituents into the venous system.

The following additional undesirable effects have occurred after using acrylic cements: temporary reduction in blood pressure, elevated serum levels of gamma-glutamyl transferase (gamma-GT) up to 10 days after surgery, thrombophlebitis, bleeding and haematoma, loosening or dislocation of the implant, superficial or deep wound infection, trochanteric bursitis, heterotopic ossification and trochanteric detachment, cardiovascular reactions such as temporary heart rhythm disorders, short-term cardiac conduction irregularities, arrhythmia, myocardial infarction and cardiac arrest, hypoxaemia, bronchospasm, pulmonary embolism, apoplexy. There have been rare reports of hypotension with anaphylaxis including anaphylactic shock associated with cardiac arrest and sudden death.

Other adverse events that may be linked to the use of PMMA bone cements are: allergic pyrexia, haematuria, dysuria, bladder

fishales, local neuropathy and vascular erosion or occlusion as well as postoperative irritation of the sciatic nerve because bone cement has been placed outside the intended area of application.

Monomer vapours can irritate the respiratory tract and eyes and may damage organs.

Warnings and precautions

During and immediately after application of the bone cement/implantation, blood pressure, pulse and respiration must be monitored carefully and appropriate measures taken if significant changes occur. If the patient develops pulmonary or cardiovascular symptoms, appropriate monitoring of blood loss is required. In the event of acute respiratory failure, anaesthesiological measures should be initiated immediately.

Polymerisation (hardening) of BonOs® R NF is an exothermic reaction. The heat given off during this reaction may damage bone or other tissue in the region of the implant.

The training and experience of the surgeon are very important when handling acrylic bone cements. Instructions for handling and mixing the cement and preparing the implant site must be followed carefully. Before using BonOs® R NF, the surgeon must be thoroughly familiar with its properties and handling characteristics. As the handling and setting characteristics of BonOs® R NF depend on the temperature and mixing technique, they are best determined by the surgeon on the basis of actual experience. For this reason, the surgeon is strongly recommended to carry out a trial run of the entire mixing, handling and setting process before performing a surgical procedure with BonOs® R NF.

Inadequate fixation and unexpected postoperative events impair the interface between the cement and bone. This can cause micro-movements, which may lead to the formation of a layer of fibrous tissue and premature implant failure. Premature loosening of the implant in the cement mantle is also possible. Therefore, long-term regular follow-up examinations are recommended for all patients.

Methyl methacrylate is a slightly volatile, flammable liquid. The vapours produced during the mixing process can irritate the respiratory tract and eyes and cause general malaise and headache. Such symptoms can be reduced with adequate ventilation or by using closed mixing systems. The monomer (methyl methacrylate) is lipid-soluble. Direct skin contact with the liquid monomer should be avoided as far as possible, as allergic reactions (contact dermatitis) cannot be ruled out. It is therefore advisable to wear an additional pair of polyethylene (PE) gloves under normal surgical gloves when handling the cement. The following materials have also proved to be suitable for protective gloves: PVP (polyethylene, ethylene vinyl alcohol, polyethylene) and viton butyl.

Contact lenses should be protected from the escaping monomer vapours.

Removal of bone cement in case of revision

For revision, a radical debridement is necessary to remove the loosened prosthesis, the cement and devitalized tissues. Reaming of the medullary canal followed by a thorough lavage using a pulse lavage system is recommended. The procedure to remove the implant is a case by case decision by the responsible surgeon depending on the condition of the implant and the patient.

Storage

BonOs® R NF should be stored away from direct sunlight. BonOs® R NF is stored at temperatures between 0 °C (32 °F) and +25 °C (77 °F). Do not use BonOs® R NF after the expiry date.

Shelf life/stability

The expiry date is printed on the outer carton, protective aluminium packaging and inner pouch. BonOs® R NF must not be used after this date. The contents of unused open or damaged packs must not be re-sterilized and should therefore be discarded. The polymer powder must not be used if it shows yellow discolouration. The liquid monomer, the ampoule itself, the inside of the blister, the cement powder and the inner PE/paper pouch are all sterile. The product must not be used if the packaging is damaged.

SINGLE-USE

BonOs® R NF must never be reused. Due to the functionality (hardening) of PMMA bone cements BonOs® R NF is only suitable for use within the specified application period. One unit of BonOs® R NF must be used for only one single patient.

Safe disposal

1. Mixed cement should be set before it is discarded with hospital waste.
2. For separate disposal of powder and liquid, please ask local waste disposal authorities.

Information

For further information, please contact your supplier or the manufacturer directly.

Symbol definitions

	Manufacturer		Consult instructions for use
	Use-by date		Caution
	Batch code		CE marking of conformity
	Catalogue number		Registered trademark
	Sterilized using aseptic processing techniques		Flammable
	Sterilized using ethylene oxide		Harmful
	Do not re-sterilize		Health Hazard
	Do not use if package is damaged		Ampoule
	Keep away from sunlight		Blister
	Temperature limit		Pouch
	Do not re-use		

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